SERVOX® Rota Tube duo

basic | voice | cuff | cuff voice

DE

EN

Verwendete Symbole / Symbols used.....35/36

Bestellinformationen / Order informations

Größe Size	SERVOX® Rota Tube duo basic		SERVOX® Rota Tube duo voice		SERVOX® Rota Tube duo cuff		SERVOX® Rota Tube duo cuff voice	
	PZN	REF	PZN	REF	PZN	REF	PZN	REF
7	10065791	29353	10065845	29357	10065928	29361	10065880	29365
8	10065816	29354	10065851	29358	10065934	29362	10065897	29366
9	10065822	29355	10065868	29359	10065940	29363	10065905	29367
10	10065839	29356	10065874	29360	10065957	29364	10065911	29368

Größe Size	Innenkanüle ungefenstert Inner cannula unfenestrated		Innenkanüle gefenstert Inner cannula fenestrated	
	PZN	REF	PZN	REF
7	16933509	29369	16933455	29373
8	16933515	29370	16933461	29374
9	16933544	29371	16933478	29375
10	16933550	29372	16933490	29376

CE 0197 MD MR

Orbisana Healthcare GmbH
Biberweg 24-26 | 53842 Troisdorf | Germany
Tel. +49 (0)2241 9322-0 | Fax +49 (0)2241 9322-277
info@orbisana.de | www.orbisana.de

40859/Version 2025-05-22

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INSTRUCTIONS FOR USE

SERVOX® Rota Tube duo basic | voice | cuff | cuff voice

Contents

	Foreword	20
1	Intended purpose	20
2	Indication	20
3	Contraindication	21
4	Side effects	21
5	Safety instructions	21
6	Scope of delivery	24
6.1	Cannula card	25
7	Description	25
7.1	SERVOX® Rota Tube duo voice / SERVOX® Rota Tube duo cuff voice (cannulas with fenestration)	25
7.2	SERVOX® Rota Tube duo cuff / SERVOX® Rota Tube duo cuff voice (cannulas with cuff)	26
7.3	Tracheal cannula connections	26
8	Tracheal cannula accessories	26
8.1	SERVOX® Cap	27
8.2	SERVOX® RW15 separating wedge	27
9	Application	27
9.1	Preparing and checking tracheal cannulas without cuff	27
9.2	Preparing and checking tracheal cannulas with cuff	28
9.2.1	Checking the cuff and safety valve	28
9.3	Using the 15 mm standard connector	28
9.4	Inserting the outer cannula without cuff	29
9.5	Inserting the outer cannula with cuff	29
9.6	Inserting the inner cannula	30
9.7	Changing the inner cannula	30
9.8	Changing the outer cannula without cuff	30
9.9	Changing the outer cannula with cuff	30
10	Cleaning and disinfection of tracheal cannulas and accessories	31
10.1	Safety instructions for cleaning and disinfection	31
10.2	Required aids and accessories	31
10.3	Cleaning instructions	32
10.4	Instructions for disinfection	32
11	Storage and transport	33
12	Period of use and disposal	33
13	Technical data	34
13.1	Parameters	34
14	Disclaimer	34

EN

40859/Version 2025-05-22

Foreword

These instructions for use are intended for patients, relatives and medical staff. The tracheal cannula is suitable for all age groups for use at home and in a clinical setting, depending on the selected size. The size, shape and type of tracheal cannulas must be specified and prescribed by the attending physician. The patient can use the product independently but must be given instructions on its use by medical staff beforehand. Tracheal cannulas and their accessories should be used in line with the descriptions in these instructions for use and in accordance with the recognised rules of medicine.

Tracheal cannulas and their accessories are assigned a batch number for traceability purposes. The batch number is indicated on the packaging next to "LOT". If you have any questions or problems with your product, please contact Orbisana Healthcare GmbH and provide the batch number.

Note

These instructions for use contain important information for safe and proper use. Keep them in a safe place for future reference.

1 Intended purpose

All SERVOX® Rota Tube duo versions are used in an artificial opening in the trachea (tracheostoma) in order to keep the tracheostoma open and to stabilise it. The SERVOX® Rota Tube duo voice and SERVOX® Rota Tube duo cuff voice can also be used for speaking in combination with a speech valve by patients with a voice option.

The SERVOX® Rota Tube duo cuff is suitable for use with ventilated tracheotomised patients. The SERVOX® Rota Tube duo cuff and SERVOX® Rota Tube duo voice can also be used for aspiration management.

Note

Tracheal cannulas and their accessories are intended for **reuse on a single patient** and must not be reused on other persons. There is a risk of infection.

2 Indication

- **SERVOX® Rota Tube duo basic:**
Use for tracheotomised, spontaneously breathing patients
- **SERVOX® Rota Tube duo voice:**
Use for tracheotomised, spontaneously breathing patients
- **SERVOX® Rota Tube duo cuff and SERVOX® Rota Tube duo cuff voice:**
Use for tracheotomised, spontaneously breathing and ventilated patients

3 Contraindication

- Tracheostoma is too small or in an unsuitable position
- Abnormal anatomy or pathology of the upper respiratory tracts
- Unusually deep trachea (e.g. in obese patients)
- Known allergies to polyurethane (PU)
- No use as a paediatric cannula
- In laryngectomised patients, the tracheal cannula must not be used with a speech valve.

4 Side effects

- Contamination of the surrounding tissue due to colonisation of the cannulas with e.g. yeast fungus (candida), bacteria, etc.
- Formation of haematomas (particularly during first use)
- Formation of granulation tissue or blocking of the fenestration in fenestrated cannulas with heavy or stringy secretions
- Risk of complications with tracheal cannulas with cuff e.g. pressure necrosis due to excessive cuff pressure or aspiration or loss of air during ventilation with a low cuff pressure
- When using fenestrated tracheal cannulas for ventilation, the fenestration may cause leaks in the ventilation system

EN

Note

Wearing a tracheal cannula may encourage an overprotective posture. It may also cause or exacerbate swallowing difficulties. If you have any questions or if complications occur, please contact your doctor.

5 Safety instructions

Safe and proper use of your tracheal cannula is important for your health. Ignoring these safety instructions can have significant effects on your health or even cause life-threatening injuries.

Attention: The word "Attention" indicates a situation that can cause damage to your product if not avoided.

In the following, warnings on the safe use of the product are indicated with the ⚠ symbol.

Danger: The word "Danger" indicates a likely situation that can result in grave injury or even death if not avoided.

Warning: The word "Warning" indicates a situation that, if not avoided, can result in grave injury or even death.

Caution: The word "Caution" indicates a situation that can result in minor or moderate injury if not avoided.

Please observe the following safety instructions:

Danger

Risk of suffocation: If the tracheal cannula is sealed, e.g. by secretions, or if aids such as artificial noses or speech valves are connected and these are sealed, e.g. by secretions, you may experience respiratory distress or suffocate.

- Check the cannula regularly! Change if necessary!
- Regularly suction secretions from your cannula!
- Regularly check the attached aids! Change if necessary!

Danger: It is not possible to exhale through a speaking valve. Laryngectomised patients without a voice prosthesis/shunt valve cannot breathe when their tracheostoma is sealed. They may experience respiratory distress or suffocate, injure themselves or die.

- Do not attach speech valves intended for tracheotomised patients when treating laryngectomised patients!
- Only use speech valves in patients who are awake and responsive!
- Never block speech valves! Remove blocked or clogged speech valves immediately!
- Never block the tracheostoma! Do not insert any objects into the cannula!

Risk of aspiration: If you insert a speaking valve directly into your tracheostoma, you might aspirate the speaking valve. You may injure yourself, experience respiratory distress, suffocate or die.

- Only attach speaking valves to the 15 mm standard connector!

Risk of injury with tracheal cannulas with cuff: If you place a speaking valve on a blocked tracheal cannula, you will not be able to exhale. You may experience respiratory distress, suffocate or die.

- In the case of a narrowed trachea, only attach speaking valves to fenestrated tracheal cannulas!
- Unblock tracheal cannulas before attaching a speaking valve!

Risk of injury/extubation/suffocation: When using a ventilator, the ventilation tube may become detached from the cannula or exert too much tension on the cannula. The tracheal cannula can change its position or be pulled out of the tracheostoma. You may injure yourself, experience respiratory distress or suffocate.

- Always secure the tracheal cannula with a tube holder!
- Ensure that the outer and inner cannula are firmly connected!
- Secure the connection between the tracheal cannula and ventilation tube!
- During ventilation, ensure that no pressure is exerted on the ventilation tube!

Risk of injury with tracheal cannulas with cuff: During ventilation, oxygen can diffuse into the cuff. The cuff pressure can rise or drop.

- Regularly check the cuff pressure.

Warning

Risk of injury with tracheal cannulas with cuff: Excessive cuff pressure can cause nerve damage or damage to the trachea.

- Check the cuff pressure regularly (approx. every 2 – 3 hours).

Risk of injury: When using local painkillers/anaesthetics, the feeling for the position of the tracheal cannula or the cuff may be impaired.

- Check the cannula position regularly!

Risk of aspiration with tracheal cannulas with cuff: Secretions can collect above a blocked cuff.

- Suction any secretions before and during unblocking of the cannula!

Risk of injury/aspiration: Objects, such as cleaning brushes, that are inserted into a tracheostoma can be aspirated; you could injure yourself, get into respiratory distress or suffer an infection.

- Never insert objects into your tracheostoma!
- Only clean tracheal cannulas outside the body!

Risk of inflammation/irritation: Residues of cleaning agents and disinfectants on the tracheal cannula can lead to irritation or inflammation.

- Rinse the tracheal cannula thoroughly after cleaning and disinfection!

Risk of aspiration: Immediately after a meal, food residues can remain in the throat and mouth and be aspirated.

- Do not change/insert tracheal cannulas until approx. 30 minutes after a meal

Risk of injury: If a tracheal cannula shows signs of damage such as alterations, discolouration or damage to the surfaces or dents or creases in the cuff, you could injure yourself.

- Do not use a damaged or altered tracheal cannula!
- Replace the tracheal cannula immediately even in the event of minor damage!
- Do not modify the tracheal cannula yourself!

Risk of injury: When attaching or removing aids, such as artificial noses or speech valves, from an inserted tracheal cannula, you may injure your trachea by exerting too much pressure or tension on the tracheal cannula.

- Hold the neck flange firmly during attachment or removal!
- Adjust pressure or tension carefully!
- If applicable, use a SERVVOX® separating wedge!

Risk of respiratory distress: Inserting a tracheal cannula increases breathing resistance. Aids such as speech valves or artificial noses also increase breathing resistance.

- If respiratory distress suddenly occurs, remove all aids immediately!
- Suction secretions if necessary. If symptoms persist, consult a doctor as soon as possible!

EN

Warning

Danger of respiratory distress when using a speaking valve: If the valve flap does not open or if, for example, the speaking valve or the valve flap is blocked or clogged with secretions, you may enter respiratory distress.

- Ensure that the valve flap assembly is always on top and that the valve flap can move freely! Remove the speaking valve immediately in case of respiratory distress!
- Do not bend over for extended periods of time while wearing a speech valve!
- Do not wear the speaking valve while lying down! Never sleep with the speaking valve on!
- Clean the speaking valve at least once or twice a day!

Risk of respiratory distress/disconnection: Residual lubricant inside the cannula can make breathing difficult. Lubricant residue on the 15 mm standard connector can cause the connection of ventilators or aids to loosen.

- Carefully remove any remaining lubricant before inserting the cannula.

Risk of allergy: Lubricants for tracheal cannulas may cause allergic reactions.

- Watch out for allergic reactions, do not use lubricants in this case!
- If you have any known allergies, check the ingredients of the lubricant!

MRI safety information

Medical devices with an MR safe symbol are safe for MRI.
SERVOX® Rota Tube duo tracheal cannulas can be used during an X-ray.

Note

Please inform us or the responsible authorities if you experience a life-threatening illness or injury or permanent physical damage as a result of using our product, or require hospital treatment or have to extend hospital treatment due to using our product.

6 Scope of delivery

- **SERVOX® Rota Tube duo basic:** 1 outer cannula, 2 inner cannulas with flexible 15 mm standard connector, 1 tube holder, 1 obturator, 1 separating wedge, 1 cannula card, 2 labels and 1 instructions for use booklet.
- **SERVOX® Rota Tube duo voice:** 1 fenestrated outer cannula, 1 inner cannula non-fenestrated with flexible 15 mm standard connector, 1 inner cannula

fenestrated with flexible 15 mm standard connector, 1 tube holder, 1 obturator, 1 separating wedge, 1 cap, 1 cannula card, 2 labels and 1 instructions for use booklet.

- **SERVOX® Rota Tube duo cuff:** 1 outer cannula with low pressure cuff, 2 inner cannulas with flexible 15 mm standard connector, 1 tube holder, 1 obturator, 1 separating wedge, 1 cannula card, 2 labels and 1 instructions for use booklet.
- **SERVOX® Rota Tube duo cuff voice:** 1 fenestrated outer cannula with low pressure cuff, 1 inner cannula non-fenestrated with flexible 15 mm standard connector, 1 inner cannula fenestrated with flexible 15 mm standard connector, 1 tube holder, 1 obturator, 1 separating wedge, 1 cap 1 cannula card, 2 labels and 1 instructions for use booklet.

6.1 Cannula card

All cannulas are supplied with a cannula card for entering patient data and the name of the attending physician. One of the enclosed labels with the information and batch number of the cannula should be attached to the back. Always take the cannula card with you and have it to hand so that the necessary information about the cannula and medical treatment is available on request or in the event of an accident.

7 Description

The cannula is inserted into the windpipe (trachea). The neck flange is located at the top end of the tube of the cannula with 2 loops on the edge of the flange to secure a holder. The cannula size, length and outer dimensions are indicated on the neck flange. The labelling is green for non-fenestrated tracheal cannulas and blue for fenestrated cannulas.

The enclosed insertion aid (obturator) facilitates the insertion of the cannula and can be used to stabilise the cannula during storage or transport. The obturator must be removed from the cannula immediately after insertion. All SERVOX® Rota Tube duo variants are sterile on delivery and must not be resterilised.

The inner cannula is inserted into the outer cannula and is secured in the neck flange hooks to prevent it from slipping out. The standard connector is cloudy and transparent for non-fenestrated inner cannulas and blue for fenestrated inner cannulas. The enclosed insertion aid (obturator) facilitates the insertion of the cannula and can be used to stabilise the cannula during storage or transport. The obturator must be removed from the cannula immediately after insertion. All SERVOX® Rota Tube duo variants are sterile on delivery and must not be resterilised.

7.1 SERVOX® Rota Tube duo voice / SERVOX® Rota Tube duo cuff voice (cannulas with fenestration)

For patients with an intact larynx, the outer cannula is fenestrated at the top of the bend and the inner cannula has a large opening (fenestration) at the same point. If the cannulas are positioned correctly, the fenestration of the two can-

nulas should lie on top of each other. By deliberately closing the cannula for a short time, the exhaled air is channelled through the fenestration via the larynx and through the vocal folds, allowing it to be used for speaking. The enclosed non-fenestrated inner cannula can be used to close the fenestration of the outer cannula, e.g. for ventilation.

7.2 SERVOX® Rota Tube duo cuff / SERVOX® Rota Tube duo cuff voice (cannulas with cuff)

The tube of the outer cannula is surrounded by an inflatable (blockable) low-pressure cuff. A pilot balloon (cuff pilot) is connected to the cuff via a supply tube. A cuff pressure gauge or a syringe with a Luer connection can be connected to the valve on the cuff pilot for filling or releasing air.

To block (inflate) the cuff, it is filled with air through this valve. To unblock the cuff, the air is released. The cuff pressure is also checked using the cuff pilot. A cuff pressure gauge must always be used to check the cuff pressure.

A blocked cuff prevents or limits secretions from the upper airways from entering the lungs (aspiration) and keeps air from flowing on the outside between the tracheal wall and the cannula when breathing; blockable tracheal cannulas must therefore always be unblocked before a speaking valve or cap is attached.

The doctor determines the cuff pressure. We recommend a cuff pressure between 15 mmHg (20 cm H₂O) and 22 mmHg (30 cm H₂O). Since movements or coughing can change the cuff pressure, it must be regularly checked, around every 2 – 3 hours, or more frequently as needed. For use in patients with a voice option, the SERVOX® Rota Tube cuff voice is fenestrated above the cuff, in the upper side of the bend. When the tracheal cannula is used together with a speaking valve for speaking, the cuff must be unblocked.

7.3 Tracheal cannula connections

The 15 mm standard connector on the upper end of the tracheal cannula can be used to attach medical aids, such as speech valves or artificial noses, or to connect a ventilator.

8 Tracheal cannula accessories

Accessories designed for use with the SERVOX® Rota Tube, such as cleaning and disinfection accessories, care sets, tracheal compresses and artificial noses, are available from Orbisana Healthcare GmbH.

Designation	CPN	Order no.
SERVOX® Tracheo Tape M duo	11600080	35207
SERVOX® RW15	–	29100
SERVOX® cap C15	02066316	29842

8.1 Cap

The cap is inserted to wean the patient off the tracheal cannula and may only be used in consultation with the doctor. Because the cap blocks breathing through the tracheostoma, the cap may only be fitted to tracheostomised patients in combination with fenestrated tracheal cannulas. When used with a cannula with cuff, the cuff must be unblocked before attaching the cap.

8.2 SERVOX® RW15 separating wedge

The SERVOX® RW15 separating wedge makes it easier to remove aids from the 15 mm standard connector. It may only be used if clean and undamaged and can be cleaned and disinfected together with the tracheal cannula.

- Carefully place the separating wedge directly behind the aid.
- Slowly push forwards onto the standard connector with a slight levering action.

EN

9 Application

Warning

- Due to the risk of aspiration, wait at least 30 minutes after eating before inserting a new tracheal cannula.
- If necessary, suction secretions before and during cannula changes.
- Always unblock the cuff before attaching a cap or a speech valve! Danger of suffocation!
- Always unblock the cuff before changing the position of the tracheal cannula or removing the tracheal cannula! Risk of injury!
- Never bend or pull on the cuff pilot or the supply tube of the cuff!
- Regularly check the cuff pressure! We recommend a cuff pressure between 15 mmHg (20 cm H₂O) and 22 mmHg (30 cm H₂O).

Note

- Wash and disinfect hands before touching the tracheostoma or tracheal cannula.
- Practise detaching and attaching various aids before inserting the cannula to facilitate subsequent handling after the tracheal cannula has been inserted.
- Relatives or medical staff who take care of the patient should wear gloves and a face mask and document all measures taken.

9.1 Preparing and checking tracheal cannulas without cuff

- Check sterile packaging and contents. Packaging and contents must not be damaged. Do not use if damaged!

- Check the cannula size. (information on the neck flange). Do not use cannulas of the wrong size!
- Check whether the obturator is easy to remove.

Improve the glide properties:

- Apply a suitable lubricant to the outer cannula or moisten the cannula with water. Carefully remove any lubricant residue inside the cannula, on the hook for engaging the inner cannula or on the connector.
- Place the outer and inner cannula next to each other within reach on a clean surface.

9.2 Preparing and checking tracheal cannulas with cuff

- Check sterile packaging and contents. Packaging and contents must not be damaged. Do not use if damaged!
- Check the cannula size. (information on the neck flange). Do not use cannulas of the wrong size!
- Check whether the obturator is easy to remove.

9.2.1 Checking the cuff and safety valve

The cuff and cuff pilot must be checked with a cuff pressure gauge to ensure they are sealed and intact every time immediately before they are used. Cannulas with a leaky cuff must not be used.

Signs of a cuff leak include, for example:

- Damage to the cuff (holes, tears, etc.).
- Air escapes noticeably from the cuff.
- Cuff cannot be completely filled.
- Water in the supply tube, cuff or control balloon (after cleaning and disinfection!).

Then:

Improve the glide properties:

- Apply a suitable lubricant to the outer cannula or moisten the cannula with water. Carefully remove any lubricant residue inside the cannula, on the hook for engaging the inner cannula or on the connector.
- Place the outer and inner cannula next to each other within reach on a clean surface.

9.3 Using the 15 mm standard connector

Attach the aid after inserting the tracheal cannula. Take care to avoid injuring the trachea.

Fitting:

- Hold the neck flange securely and carefully push the aid onto the 15 mm standard connector. The aid must sit securely enough but still be easily removable! If it is not in the correct position, it may fall out.

Removal:

- Hold the neck flange securely and carefully pull the aid off the connector, using a separating wedge if necessary.

9.4 Inserting the outer cannula without cuff

- Clean the tracheostoma.
- Put the patient in the supine position with the upper body elevated and the head slightly tilted back.
- Insert the outer cannula, as practiced, into the tracheostoma. An urge to cough may result.
- Remove the obturator immediately!
- Check that the cannula is positioned correctly. The tracheal cannula should be in the tracheostoma with as little tension as possible.

Secure the tracheal cannula:

- Attach the tube holder to the eyelets in the neck flange and pull tight. The band must be tight enough to ensure it does not come loose but not so tight that it constricts the throat! A finger should fit between the neck and the tube holder. Adjust and check the fit.
- Wait a few minutes in case any complications such as breathing problems, the urge to cough, vomiting etc. occur. The throat may briefly swell up when changing the cannula.

9.5 Inserting the outer cannula with cuff

- Clean the tracheostoma.
- Put the patient in the supine position with the upper body elevated and the head slightly tilted back.
- Insert the outer cannula with the completely unblocked cuff, as practiced, into the tracheostoma. An urge to cough may result.
- Remove the obturator immediately!
- The tracheal cannula should be in the tracheostoma with as little tension as possible.
- Block the cuff until the trachea is adequately sealed. While doing so, check the cuff pressure using the cuff pressure gauge. We recommend a cuff pressure between 15 mmHg (20 cm H₂O) and 22 mmHg (30 cm H₂O).
- Do not change the cannula position after blocking the cuff! Risk of injury!

Secure the tracheal cannula:

- Attach the tube holder to the eyelets in the neck flange and pull tight. The band must be tight enough to ensure it does not come loose but not so tight that it constricts the throat! A finger should fit between the neck and the tube holder. Adjust and check the fit.
- Wait a few minutes in case any complications such as breathing problems, the urge to cough, vomiting etc. occur. The throat may briefly swell up when changing the cannula.
- Check the cuff pressure regularly (approx. every 2 – 3 hours).

9.6 Inserting the inner cannula

Hold the outer cannula securely by the neck flange and carefully insert the inner cannula and click into the holders. Check that it is correctly engaged.

9.7 Changing the tracheal cannula

Change and clean/disinfect the inner cannula at the interval specified by the doctor or at least 1 – 2 times a day. If secretions are deposited inside the inner cannula and cannot be removed by coughing or suction, the inner cannula must be changed and cleaned. If necessary, suction any secretions before changing the cannula.

- Prepare a replacement inner cannula.
- Hold the outer cannula securely by the neck flange, carefully release the inner cannula from the hooks and pull out.
- Suction any secretions if necessary.
- Carefully insert and secure the new inner cannula.
- Clean and disinfect the used inner cannula immediately.

9.8 Changing the outer cannula without cuff

Change the tracheal cannula at the interval specified by the doctor. Always suction secretions before changing the cannula. Clean and disinfect the tracheal cannula every 7 days (more frequently if necessary). The period of use of the tracheal cannula is 29 days. During this time, it may be cleaned and disinfected a maximum of 8 times.

- Prepare a replacement cannula and have it ready.
- Open the tube holder and slowly and carefully remove the tracheal cannula from the tracheostoma. This may cause an urge to cough.
- Clean the tracheostoma and allow to air dry.
- Insert the new outer cannula.
- Clean and disinfect the used tracheal cannula immediately.

9.9 Changing the outer cannula with cuff

Change the tracheal cannula at the interval specified by the doctor. Always suction secretions before changing the cannula. Clean and disinfect the tracheal cannula every 7 days (more frequently if necessary). The period of use of the tracheal cannula is 29 days. During this time, it may be cleaned and disinfected a maximum of 8 times.

- Prepare a replacement outer cannula and have it ready.
- Suction secretions before and during unblocking.
- Completely unblock the cuff!
- Open the tube holder and slowly and carefully remove the outer cannula from the tracheostoma. This may cause an urge to cough.
- Clean the tracheostoma and allow to air dry.
- Insert the new outer cannula.
- Clean and disinfect the used outer cannula promptly.

10 Cleaning and disinfection of tracheal cannulas and accessories

Note

It is often not possible to clean and disinfect individual tracheal cannulas for reuse on a single patient in a clinical setting or this may contradict internal hospital regulations. When used in clinical settings, it is therefore recommended to use the tracheal cannula as a single-use product and to dispose of it after use. In a domestic environment, the period of use is 29 days. During this time, the cannula can be cleaned and disinfected up to 8 times as described here.

EN

10.1 Safety instructions for cleaning and disinfection

Warning

- Only clean and disinfect tracheal cannulas and accessories outside the body!
- Always rinse tracheal cannulas and accessories thoroughly after cleaning and disinfection!
- Never insert objects into a tracheostoma!

Attention

- Only use with suitable and recommended cleaning agents or disinfectants!
- Do not use any aggressive household cleaners or highly concentrated alcohol!
- Only use suitable cleaning brushes! Do not use sharp objects such as scrapers, knives or fingernails, as these can damage the tracheal cannula!
- Clean tracheal cannulas and accessories before disinfection!
- **Never brush off the cuff!**
- Never immerse the cuff pilot in cleaning/disinfectant solution!
- Do not resterilise tracheal cannulas and accessories!
- Only recondition the tracheal cannula at a washbasin with a hot water connection! Thoroughly clean and disinfect the washbasin before reconditioning the tracheal cannula!

10.2 Required aids and accessories

- Appropriately sized cannula cleaning brush. The cleaning brush should be flexible, with soft bristles and a rounded tip in order to prevent damage to the cannula.

- Suitable plastic cleaning tub with a sieve insert (e. g. SERVOX® cleaning tub with sieve insert).
- Suitable (cannula) cleaning agent for alkaline-stable instruments and devices in all medical areas and laboratories (biodegradable surfactants: non-ionic surfactants < 5 %, anionic surfactants 5 – 15 %, phosphate > 30 %).
- Cannula disinfectant for plastic tracheal cannulas.
- Gloves.
- Access to warm, clean, and running water.

10.3 Cleaning instructions

- Wash and disinfect hands thoroughly.
- Remove and clean or dispose of attached aids (observe instructions for use).
- Inflate the cuff before cleaning.
- To remove larger impurities, carefully flush the tracheal cannulas and the accessories with lukewarm running water.
- Remove firmly attached secretions with the cannula cleaning brush.

Inside:

- Bend the cleaning brush into the shape of the cannula, insert from the tip of the cannula and carefully move up and down and turn inside with slight pressure for at least 30 seconds. Each time before removing the brush, rinse the brush head under running water.
- Brush the outside of the cannulas for a further 30 seconds. Avoid strong pressure – do not brush off the cuff!
- Then soak the tracheal cannulas and accessories in cleaning solution.
- **The cuff pilot must not be immersed in the cleaning solution!**
- **Prepare the cleaning solution with the cannula cleaning powder** in the cleaning tub according to the manufacturer's instructions. The tracheal cannulas and accessories must be completely covered during the exposure time.

Avoid skin contact with the cleaning solution:

- Wear gloves.
- Place the tracheal cannulas and accessories in the sieve insert of the cleaning tub and immerse without the cuff; do not leave any air inside the cannula. All parts must be completely immersed and soaked for at least 10 minutes (observe manufacturer instructions).
- Then carefully lift out the sieve insert and remove all parts.
- Clean the tracheal cannulas again with the cleaning brush.
- Then rinse all parts thoroughly several times under clear running water.
- Then dry all parts with a lint-free cloth or leave to dry completely at room temperature.
- Dispose of the cleaning solution in accordance with the manufacturer's instructions.

10.4 Instructions for disinfection

- Wash and disinfect hands thoroughly.
- Inflate the cuff before disinfection.

- Soak the tracheal cannula and accessories in disinfectant solution.
- The cuff pilot must not be immersed in the disinfectant solution!

Disinfectant solution:

- Prepare a 2 % disinfectant solution with the cannula disinfectant in the cleaning tub according to the manufacturer's instructions. The tracheal cannula and accessories must be completely covered during the exposure time.
- Avoid skin contact with the disinfectant solution:
- Wear gloves (observe manufacturer's instructions).
- Place the previously cleaned tracheal cannulas and accessories in the sieve insert of the cleaning tub and immerse without the cuff; do not leave any air inside the cannulas. All parts must be completely immersed and soaked for at least 30 minutes (observe manufacturer instructions).
- Then carefully lift out the sieve insert and remove all parts.
- Then rinse all parts thoroughly several times under clear running water.
- Then dry all parts with a lint-free cloth or leave to dry completely at room temperature.
- Dispose of the disinfectant solution in accordance with the manufacturer's instructions.

EN

Note

The suitability of the manual cleaning and disinfection process described here was tested by an independent, officially accredited test laboratory using the cleaning powder SERVOX® cannula cleaning powder (concentration according to manufacturer's instructions, Dr. Schumacher GmbH) and the cannula disinfectant PRIMEDICLEAN (concentration according to manufacturer's instructions, Primed Halberstadt Medizintechnik GmbH).

11 Storage and transport

Store or transport cleaned and disinfected tracheal cannulas and accessories in a clean and dry environment, ideally in a lockable container free from dust and protected from light. Protect from direct sunlight and store or transport at temperatures between 5 °C und 49 °C. For tracheal cannulas with cuff: Do not block the cuff during storage or transport.

12 Period of use and disposal

In intact sterile packaging and in compliance with the storage conditions, tracheal cannulas and accessories have a maximum shelf life of 3 years. Tracheal cannulas and accessories can be used multiple times on one patient for a maximum of 29 days after opening the packaging, provided they have been properly cleaned and disinfected a maximum of 8 times and the instructions for use have been followed.

Damaged tracheal cannulas or accessories must not be used, they must be disposed of. Signs include alterations, discolouration or damage to the surfaces of accessories, cannula or cuff as well as dents or creases in the cuff. All tracheal cannulas and their accessories must be disposed of in line with national regulations for waste products. If no special measures are necessary due to illness (e.g. MRSA infection), tracheal cannulas and their accessories can be disposed of with household waste.

13 Technical data

All versions of the SERVOX® Rota Tube duo have a bending angle of 90°. The outer cannulas are made of polyurethane, and the inner cannulas are made of polypropylene and contain no metallic materials. There are no restrictions on the use of magnetic resonance imaging (MRI).

The SERVOX® Rota Tube duo cuff and SERVOX® Rota Tube duo cuff voice variants are fitted with radiopaque strip on the supply tube to check the position.

13.1 Parameters


Diagram of the tracheal cannula

Size	Outer diameter of outer cannula OD* (in mm)	Inner diameter of inner cannula ID ** (in mm)	Length L*** (in mm)	Angle θ (in °)	Cuff diameter **** (in mm)
7	10.6	7	76	90	23
8	11.7	8	78	90	26
9	12.8	9	80	90	30
10	13.8	10	82	90	33

* Outer diameter measured at the widest cross-section, ** Inner diameter measured at the smallest inner lumen, *** Length measured over the centre arch, **** Only for version with cuff

14 Disclaimer

No repairs or changes may be made to the product by persons not authorised by Orbisana Healthcare GmbH. In the event of damage, the manufacturer does not assume liability in case of non-compliance with the instructions for use or the safety instructions or manipulation of the product.

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