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41729/Version 2022-03-22

# Instructions for use SERVOX® Tracheo Nose F

Artificial nose with pressure relief valve

REF 25732

## 1. General description

The SERVOX® Tracheo Nose F is an artificial nose with a 15 mm-standard connector. It can be attached to endotracheal tubes and tracheostomy tubes with a 15 mm-standard connector. The foam material inside the artificial nose moistens, warms and filters the inspired air. It has a 5 mm oxygen port on its side to which an oxygen safety tube can be connected. The PVC pressure relief valve (transparent cap) prevents the artificial nose from getting detached from the tube during coughing and prevents an excessive increase in pressure in the respiratory tract.

## 2. Indication

The SERVOX® Tracheo Nose F is intended for spontaneously breathing tracheotomized or laryngectomized patients with a tracheostomy tube. It can also be used for spontaneously breathing patients who breathe through an endotracheal tube. The SERVOX® Tracheo Nose F is suitable for patients with a tidal volume of 250 – 1000 ml. The SERVOX® Tracheo Nose F can be used in clinics or home care. The SERVOX® Tracheo Nose F is latex-free and does not contain any metal materials therefore no restrictions on the use of magnetic resonance imaging (MRI) apply. It is electrically non-conducting.

The SERVOX® Tracheo Nose F is supplied non-sterile and is intended for single use only.

## 3. Contraindications

The SERVOX® Tracheo Nose F is to be used only within the recommended tidal volume. If the tidal volume is too low, the dead space increases which can lead to an increase in CO<sub>2</sub> concentration. A too high tidal volume results in insufficient moistening of the air. The use of the SERVOX® Tracheo Nose F is not recommended for the following patients:

- Patients with severe exsiccosis (dehydration)
- Patients with excessive secretion production in the respiratory tract
- Patients using a speech valve

## 4. Side effects

No side effects are known, in cases where account has been taken of the indication, contraindications and instructions for use.

## 5. Application

When attaching the SERVOX® Tracheo Nose F to the 15 mm standard connector of a tracheostomy tube or an endotracheal tube, it should be ensured that the device is firmly fitted to the 15 mm-standard connector.

An oxygen safety tube can be connected to the oxygen port on the side. It should be noted that the moistening and warming effect begins only after a few breaths. The SERVOX® Tracheo Nose F should be replaced after 24 hours at the latest. A more frequent change may be indicated in cases of severe mucous production or soiling of the filter medium.

## 6. Warnings

- If there is an increase in airway resistance caused for example by amplified accumulations of secretions in the SERVOX® Tracheo Nose F, the artificial nose has to be removed at once, disposed of immediately and replaced by a new SERVOX® Tracheo Nose F.
- Used artificial noses are to be disposed of immediately.

Failure to comply with these instructions can lead to severe, even to life-threatening damage to the patient's health. The manufacturer undertakes no liability in the event of injury arising from failure to comply with these instructions and tampering with the product.

## 7. Safety

A batch number is assigned to every product, which is shown on the packaging and label in order to ensure clear traceability.

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### Symbol explanation

	CE-mark with number of notified body		Catalogue number
	Consult instructions for use		Batch code
	Do not re-use		Pharmacy Central Number
	Keep away from sunlight		The green dot
	Keep dry		Medical device
	Do not use if package is damaged		Temperature limit storage and operation +5°C - +25°C
	Date of manufacture		Non-sterile
	Manufacturer		Humidity limitation storage and operation ≤ 65% RH

Servona GmbH

Biberweg 24-26 | D-53842 Troisdorf  
Tel. 0 22 41.93 22-0  
Fax 0 22 41.93 22-277  
info@servona.de | www.servona.de

CE 0197