SERVOX®

HME Modul

comfort | comfort skin

HME Modul+
HF | RF

DE

EN

ES

FR

IT

NL

HR

LT

PL

44926/Version 2025-04-28

Orbisana Healthcare GmbH

Biberweg 24-26
53842 Troisdorf | Germany
Tel. +49(0)2241 9322-0
Fax. +49(0)2241 9322-277
info@orbisana.de
www.orbisana.de

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SERVOX® HME Modul

comfort | comfort skin

SERVOX® HME Modul+ HF | RF

Verwendete Symbole	47
--------------------------	----

Inhalt

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Contents

1 Foreword	7
1.1 Intended purpose	8
1.2 Indications	8
1.3 Contraindications	8
2 Safe use	8
3 Description	9
4 Application	10
5 Technical data	10
6 Disclaimer	11

1 Foreword

Name	colour	PU	Order no.
SERVOX® HME Modul+ HF	white	30	38317
SERVOX® HME Modul+ RF	white	30	38318
SERVOX® HME Modul comfort	white	30	29454
SERVOX® HME Modul comfort skin	beige	30	30015

These instructions for use are intended for patients, relatives and healthcare professionals. The product can be used in home care and clinical settings. Patients

can use the product unaided. Use the product according to these instructions for use and according to the recognised rules of medicine.

For clear traceability, a batch code is assigned to the product. It can be found on the packaging behind the „LOT“ marking. If you have any questions or problems with the product, please contact the sales representative or Orbisana Healthcare GmbH, quoting the batch code.

Note

These instructions for use contain important information for safe and proper use. Keep it for future reference.

1.1 Intended purpose

The SERVOX® HME Modul humidifies, warms and cleans breathing air while restoring natural airway resistance as far as possible. It also enables targeted closing of the tracheostoma and thus phonation at the press of a button.

Note

The SERVOX® HME Modul is a nonsterile, disposable medical product and must not be reused, cleaned or reprocessed.

1.2 Indications

- For spontaneously breathing tracheotomised patients with largely preserved glottis (vocal folds) who can perform phonation by closing the tube with their fingers.
- For spontaneously breathing laryngectomised patients, phonation use is only possible when using a shunt valve (vocal prosthetic).

1.3 Contraindications

- Patients without spontaneous breathing
- In case of strong secretion production
- In case of secretion and water accumulation in the lungs and airways (left heart failure, pulmonary oedema),
- In case of strong air leakage (exhaled air does not flow through the HME).
- In case of insufficient tidal volume, as the additional dead space could cause impaired gas exchange

2 Safe use

Safe and proper handling is important for your health! Ignoring the safety instructions can have significant effects on your health or even cause life-threatening situations.

Warning: The word “Warning” indicates a situation that, if not avoided, can result in grave injury or even death.

Please heed the following safety instructions:

Warning

Risk of infection: Using a soiled SERVOX® HME Modul puts you at risk of infection.

- Exchange the SERVOX® HME Modul within 24 hours!
- Exchange the SERVOX® HME Modul more frequently in case of heavy mucus secretion or soiling of the filter medium!
- Do not use the SERVOX® HME Modul if the SERVOX® HME Modul or its packaging is damaged!
- Do not use a SERVOX® HME Modul that has fallen out of the holder or which was accidentally ejected (e.g. due to severe coughing)! Replace it with a new SERVOX® HME Modul!
- Do not rinse the SERVOX® HME Modul with water!
- Dispose of used SERVOX® HME Modul units immediately after use!

Risk of suffocation / risk of respiratory distress: If your tracheostoma is blocked accidentally or continuously, you are at risk of respiratory distress or suffocation.

Laryngectomised patients without shunt valves cannot breathe when their tracheostoma is closed.

- Never place the SERVOX® HME Modul directly in your tracheostoma!
- Do not make any changes to the SERVOX® HME Modul!
- Replace the SERVOX® HME Modul immediately in case of clogged filter or increased airway resistance e.g. due to increased accumulation of secretion!
- Make sure not to block the front or the bottom of the SERVOX® HME Modul with clothing, etc.!
- Do not use any atomisers or humidifiers and do not add any water to the SERVOX® HME Modul, as this can increase the breathing resistance and block respiration!
- Do not use the SERVOX® HME Modul in water, the shower or the bathtub!

EN

3 Description

The SERVOX® HME Modul is a heat and moisture exchanger (HME). Moisture and heat of the expired air are preserved in a filter medium upon exhale and then returned to the breathing air on inhale. The humidifying effect sets in after just a few breaths. Moistening and heating of the breathing air minimises complications that can occur when patients lose their capacity to breathe through the upper airways. Using the SERVOX® HME Modul increases the airway resistance, similar to the physiological resistance of the upper airways, supporting natural lung function.

The SERVOX® HME Modul consists of a plastic housing with a filter medium. On the patient side, there is a 22 mm connector **A** for inserting the SERVOX® HME Modul into auxiliaries with a 22 mm holder, such as a tracheal cannula, a

44926/Version 2025-04-28

base plate or a stoma button. The front of the SERVVOX® HME Modul features a button [B] labelled "SERVOX®" that is suspended by a spring element.

Figure 1: SERVVOX® HME Modul diagram*

*Note: The original product may appear different from diagram.

4 Application

Checking the SERVVOX® HME Modul: Check the packaging and product for damage.

Inserting the SERVVOX® HME Modul: With clean hands, lightly press the SERVVOX® HME Modul with the patient-side 22 mm connector into the 22 mm holder of the tracheal cannula, base plate or stoma button. The SERVVOX® HME Modul must sit tight in the 22 mm holding ring! If the module is not inserted properly, it may fall out.

Using the SERVVOX® HME Modul: Pressing the button on the front of the SERVVOX® HME Modul with the finger compresses the spring element, closing the SERVVOX® HME Modul. The airway through the tracheostoma is blocked. In tracheotomised patients, the expired air is lead through the natural airways; in laryngectomised patients, it is guided through the shunt valve into the oesophagus. In this way, expired air can be used for phonation. The tracheostoma can also be closed to increase pressure in the airways, e.g. to cough.

Removing the SERVVOX® HME Modul: With clean hands, hold the 22 mm holder of the tracheal cannula, base plate or stoma button and carefully pull the SERVVOX® HME Modul out of the 22 mm holder. Dispose of the SERVVOX® HME Modul immediately after use.

5 Technical data

The SERVOX® HME Modul is made of latex-free plastic, meaning there are no restrictions in performing a magnetic resonance imaging (MRI) examination on a patient.

The SERVOX® HME Modul is a nonsterile, disposable medical product and must not be reused, cleaned or reprocessed. The SERVOX® HME Modul can be combined with auxiliaries via a 22 mm holder.

	comfort/ comfort skin	+RF (Regular Flow)	+HF (High Flow)
Tidal volume:	50 – 1.000 ml		
Inner volume/dead space:	5,7 ml		
Max. period of use:	24 hours		
Connection:	22 mm connecting piece for use in 22 mm holder		
Weight:	2,7 g	2,6 g	

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Keep dry and protect from sunlight. From the date of manufacture, the SERVOX HME Modul has a shelf life of 3 years if the package is unopened and stored correctly.

6 Disclaimer

In the event of damage, the manufacturer does not assume liability in case of non-compliance with the instruction manual or the safety instructions or manipulation of the product.

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